


GNE-3511

Chemical Properties

CAS No. : 1496581-76-0
 Formula: C₂₃H₂₆F₂N₆O
 Molecular Weight: 440.49
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	GNE-3511 is a potent and selective dual leucine zipper kinase inhibitors.
Targets(IC50)	MAPK,DNA Alkylation

Solubility Information

Solubility	DMSO: 6.25 mg/mL (14.19 mM),Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
In vivo Formulation	10% DMSO+40% PEG300+5% Tween 80+45% Saline: 1 mg/mL (2.27 mM),Sonication is recommended. <i>Please add the solvents sequentially, clarifying the solution as much as possible before adding the next one. Dissolve by heating and/or sonication if necessary. Working solution is recommended to be prepared and used immediately. The formulation provided above is for reference purposes only. In vivo formulations may vary and should be modified based on specific experimental conditions.</i>

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2702 mL	11.351 mL	22.702 mL
5 mM	0.454 mL	2.2702 mL	4.5404 mL
10 mM	0.227 mL	1.1351 mL	2.2702 mL
50 mM	0.0454 mL	0.227 mL	0.454 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Patel S, et al. J Med Chem. 2015 Jan 8;58(1):401-18.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481