


2'-Fluoro-5MeC (Bz)-3'-phosphoramidite

Chemical Properties

CAS No. : 182495-83-6
 Formula: C47H53FN5O8P
 Molecular Weight: 865.92
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	2'-Fluoro-5MeC (Bz)-3'-phosphoramidite is a Nucleoside Phosphoramidite.
Targets(IC50)	Nucleoside Antimetabolite/Analog,Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.1548 mL	5.7742 mL	11.5484 mL
5 mM	0.231 mL	1.1548 mL	2.3097 mL
10 mM	0.1155 mL	0.5774 mL	1.1548 mL
50 mM	0.0231 mL	0.1155 mL	0.231 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481