


γ -Secretase modulator 11

Chemical Properties

CAS No. :	2434630-29-0
Formula:	C ₂₈ H ₂₂ F ₂ N ₄ O ₂
Molecular Weight:	484.5
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	5-{8-[(3,4'-Difluoro[1,1'-biphenyl]-4-yl)methoxy]-2-methylimidazo[1,2-a]pyridin-3-yl}-N-methylpyridine-2-carboxamide (1o) exhibited high in vitro potency and brain exposure, significantly reduced brain A β 42 levels, showed no inhibitory effect on cytochrome p450 enzymes, and demonstrated excellent anti-cognitive deficits in AD model mice.
Targets(IC50)	Others,Gamma-secretase

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.064 mL	10.3199 mL	20.6398 mL
5 mM	0.4128 mL	2.064 mL	4.128 mL
10 mM	0.2064 mL	1.032 mL	2.064 mL
50 mM	0.0413 mL	0.2064 mL	0.4128 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481