


2'-Deoxy-isoguanosine

Chemical Properties

CAS No. :	106449-56-3
Formula:	C10H13N5O4
Molecular Weight:	267.24
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	2'-Deoxy-isoguanosine is a Nucleoside Derivative - Iso-guanosine Derivative.
Targets(IC50)	Nucleoside Antimetabolite/Analog,Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.742 mL	18.7098 mL	37.4195 mL
5 mM	0.7484 mL	3.742 mL	7.4839 mL
10 mM	0.3742 mL	1.871 mL	3.742 mL
50 mM	0.0748 mL	0.3742 mL	0.7484 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481