

Tecomin

Chemical Properties


CAS No. :

Formula: C₁₅H₂₀O₉

Molecular Weight: 344.316

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

| | |
|---------------|---|
| Description | Tecomin is a useful organic compound for research related to life sciences and the catalog number is T125081. |
| Targets(IC50) | Others |

Preparing Stock Solutions

| | 1mg | 5mg | 10mg |
|-------|-----------|------------|------------|
| 1 mM | 2.9043 mL | 14.5214 mL | 29.0428 mL |
| 5 mM | 0.5809 mL | 2.9043 mL | 5.8086 mL |
| 10 mM | 0.2904 mL | 1.4521 mL | 2.9043 mL |
| 50 mM | 0.0581 mL | 0.2904 mL | 0.5809 mL |

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481