

SR 19881

Chemical Properties


CAS No. : 2213490-89-0

Formula: C₁₉H₂₄N₂O₂

Molecular Weight: 312.41

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	SR 19881 is a full agonist of ERRγ (EC ₅₀ value of 0.39 μM in a binding assay and an EC ₅₀ value of 4.7 μM in a cell-based assay).
Targets(IC ₅₀)	Estrogen Receptor/ERR, Estrogen/progestogen Receptor

Solubility Information

Solubility	DMSO: 22.5 mg/mL (72.02 mM), Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
In vivo Formulation	10% DMSO+40% PEG300+5% Tween-80+45% Saline: 2 mg/mL (6.4 mM), Sonication is recommended. <i>Please add the solvents sequentially, clarifying the solution as much as possible before adding the next one. Dissolve by heating and/or sonication if necessary. Working solution is recommended to be prepared and used immediately. The formulation provided above is for reference purposes only. In vivo formulations may vary and should be modified based on specific experimental conditions.</i>

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.2009 mL	16.0046 mL	32.0092 mL
5 mM	0.6402 mL	3.2009 mL	6.4018 mL
10 mM	0.3201 mL	1.6005 mL	3.2009 mL
50 mM	0.064 mL	0.3201 mL	0.6402 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Lin H , Doebelin C , Patouret, Rémi, et al. Design, synthesis, and evaluation of simple phenol amides as ERRy agonists[J]. Bioorganic & Medicinal Chemistry Letters, 2018:S0960894X18301987.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481