

## Ganoderic acid G

## Chemical Properties

CAS No. : 98665-22-6


Formula: C<sub>30</sub>H<sub>44</sub>O<sub>8</sub>

Molecular Weight: 532.67

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Ganoderic acid G is isolated from the surface of Ganoderma lucidum tube layer.
Targets(IC50)	Others

## Solubility Information

Solubility	DMSO: 12 mg/mL (22.53 mM),Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
In vivo Formulation	10% DMSO+40% PEG300+5% Tween 80+45% Saline: 2 mg/mL (3.75 mM),Sonication is recommended. <i>Please add the solvents sequentially, clarifying the solution as much as possible before adding the next one. Dissolve by heating and/or sonication if necessary. Working solution is recommended to be prepared and used immediately. The formulation provided above is for reference purposes only. In vivo formulations may vary and should be modified based on specific experimental conditions.</i>

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.8773 mL	9.3867 mL	18.7733 mL
5 mM	0.3755 mL	1.8773 mL	3.7547 mL
10 mM	0.1877 mL	0.9387 mL	1.8773 mL
50 mM	0.0375 mL	0.1877 mL	0.3755 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Ganoderic acid G and I and ganolucidic acid A and B, new triterpenoids from *Ganoderma lucidum*.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481