


Peperomin C

Chemical Properties

CAS No. :	129341-10-2
Formula:	C ₂₄ H ₃₀ O ₈
Molecular Weight:	446.496
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Peperomin C is a useful organic compound for research related to life sciences. The catalog number is T126092 and the CAS number is 129341-10-2.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2396 mL	11.1982 mL	22.3964 mL
5 mM	0.4479 mL	2.2396 mL	4.4793 mL
10 mM	0.224 mL	1.1198 mL	2.2396 mL
50 mM	0.0448 mL	0.224 mL	0.4479 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481