


BITMC

Chemical Properties

CAS No. :	1622922-42-2
Formula:	C ₂₀ H ₁₄ BrN ₅ O ₅
Molecular Weight:	452.33
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	BITMC is a novel potent IDO1 inhibitor.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2108 mL	11.0539 mL	22.1078 mL
5 mM	0.4422 mL	2.2108 mL	4.4216 mL
10 mM	0.2211 mL	1.1054 mL	2.2108 mL
50 mM	0.0442 mL	0.2211 mL	0.4422 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481