


## Barbinervic acid

## Chemical Properties

CAS No. :	64199-78-6
Formula:	C <sub>30</sub> H <sub>48</sub> O <sub>5</sub>
Molecular Weight:	488.7
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Barbinervic acid shows anti-bacteria activity (IC <sub>50</sub> = 0.098 g x LU(-1), IC <sub>50</sub> = 0.27 g x L (-1)) against SA and ESBLs-SA respectively.
Targets(IC <sub>50</sub> )	Glutaminase,Antifection

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0462 mL	10.2312 mL	20.4625 mL
5 mM	0.4092 mL	2.0462 mL	4.0925 mL
10 mM	0.2046 mL	1.0231 mL	2.0462 mL
50 mM	0.0409 mL	0.2046 mL	0.4092 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481