


Wy 27569

Chemical Properties

CAS No. : 103417-69-2
 Formula: C₂₁H₂₂N₄O₆
 Molecular Weight: 426.42
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	Wy 27569 is a combined calcium channel blocker and thromboxane synthetase inhibitor.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3451 mL	11.7255 mL	23.4511 mL
5 mM	0.469 mL	2.3451 mL	4.6902 mL
10 mM	0.2345 mL	1.1726 mL	2.3451 mL
50 mM	0.0469 mL	0.2345 mL	0.469 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481