


Tokinolide B

Chemical Properties

CAS No. :	112966-16-2
Formula:	C ₂₄ H ₂₈ O ₄
Molecular Weight:	380.48
Storage:	Keep away from direct sunlight, Store at low temperature, Keep away from moisture Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	Tokinolide B, isolated from Ligusticum porter rhizomes, mediates anti-inflammatory effects via Nur77 binding.
Targets(IC50)	Others

Solubility Information

Solubility	DMSO: 80 mg/mL (210.26 mM), Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
In vivo Formulation	10% DMSO+40% PEG300+5% Tween 80+45% Saline: 3.3 mg/mL (8.67 mM), Sonication is recommended. <i>Please add the solvents sequentially, clarifying the solution as much as possible before adding the next one. Dissolve by heating and/or sonication if necessary. Working solution is recommended to be prepared and used immediately. The formulation provided above is for reference purposes only. In vivo formulations may vary and should be modified based on specific experimental conditions.</i>

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.6283 mL	13.1413 mL	26.2826 mL
5 mM	0.5257 mL	2.6283 mL	5.2565 mL
10 mM	0.2628 mL	1.3141 mL	2.6283 mL
50 mM	0.0526 mL	0.2628 mL	0.5257 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Phthalides and other constituents from *Ligusticum porteri*; sedative and spasmolytic activities of some natural products and derivatives. *Nat Prod Res.* 2011 Aug;25(13):1234-42.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481