


## Cyanidin 3-sophoroside chloride

### Chemical Properties

CAS No. : 18376-31-3  
 Formula: C<sub>27</sub>H<sub>31</sub>ClO<sub>16</sub>  
 Molecular Weight: 646.98  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
*Actual storage temperature shall be subject to the COA.*



### Biological Description

Description	Cyanidin-3-sophoroside has antioxidant capacity.
Targets(IC50)	Others, Tyrosinase

### Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.5456 mL	7.7282 mL	15.4564 mL
5 mM	0.3091 mL	1.5456 mL	3.0913 mL
10 mM	0.1546 mL	0.7728 mL	1.5456 mL
50 mM	0.0309 mL	0.1546 mL	0.3091 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481