


## Conantokin-T

## Chemical Properties

CAS No. :	127476-26-0
Formula:	C110H175N31O45S
Molecular Weight:	2683.81
Storage:	Keep away from moisture Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



## Biological Description

Description	Non-competitive NMDA receptor antagonist
Targets(IC50)	Others,iGluR

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.3726 mL	1.863 mL	3.726 mL
5 mM	0.0745 mL	0.3726 mL	0.7452 mL
10 mM	0.0373 mL	0.1863 mL	0.3726 mL
50 mM	0.0075 mL	0.0373 mL	0.0745 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481