


23-Hydroxymangiferonic acid

Chemical Properties

CAS No. :	232266-08-9
Formula:	C30H46O4
Molecular Weight:	470.694
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	23-Hydroxymangiferonic acid is a natural product for research related to life sciences. The catalog number is TN2815 and the CAS number is 232266-08-9.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1245 mL	10.6227 mL	21.2454 mL
5 mM	0.4249 mL	2.1245 mL	4.2491 mL
10 mM	0.2125 mL	1.0623 mL	2.1245 mL
50 mM	0.0425 mL	0.2125 mL	0.4249 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481