


23-Hydroxymangiferonic acid

Chemical Properties

| | |
|-------------------|---|
| CAS No. : | 232266-08-9 |
| Formula: | C ₃₀ H ₄₆ O ₄ |
| Molecular Weight: | 470.694 |
| Storage: | Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA. |



Biological Description

| | |
|---------------|---|
| Description | 23-Hydroxymangiferonic acid is a natural product for research related to life sciences. The catalog number is TN2815 and the CAS number is 232266-08-9. |
| Targets(IC50) | Others |

Preparing Stock Solutions

| | 1mg | 5mg | 10mg |
|-------|-----------|------------|------------|
| 1 mM | 2.1245 mL | 10.6227 mL | 21.2454 mL |
| 5 mM | 0.4249 mL | 2.1245 mL | 4.2491 mL |
| 10 mM | 0.2125 mL | 1.0623 mL | 2.1245 mL |
| 50 mM | 0.0425 mL | 0.2125 mL | 0.4249 mL |

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481