


## (S)-3-Phenoxypropane-1,2-diol

## Chemical Properties

CAS No. :	139165-57-4
Formula:	C <sub>9</sub> H <sub>12</sub> O <sub>3</sub>
Molecular Weight:	168.19
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	(S)-3-Phenoxypropane-1,2-diol, with CAS No. 139165-57-4, is a fragment molecule that serves as an important scaffold for molecular linking, expansion, and modification. (S)-3-Phenoxypropane-1,2-diol provides a structural basis and research tool for the design and screening of novel drug candidates, and is commonly used in drug discovery, drug synthesis, and related research.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.9457 mL	29.7283 mL	59.4566 mL
5 mM	1.1891 mL	5.9457 mL	11.8913 mL
10 mM	0.5946 mL	2.9728 mL	5.9457 mL
50 mM	0.1189 mL	0.5946 mL	1.1891 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481