

Decanal, 2-methylene-

Chemical Properties


CAS No. : 22418-65-1

Formula: C₁₁H₂₀O

Molecular Weight: 168.28

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Decanal, 2-methylene- is a biochemical.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.9425 mL	29.7124 mL	59.4248 mL
5 mM	1.1885 mL	5.9425 mL	11.885 mL
10 mM	0.5942 mL	2.9712 mL	5.9425 mL
50 mM	0.1188 mL	0.5942 mL	1.1885 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Miles WH, Dethoff EA, Tuson HH, Ulas G. Kishner's reduction of 2-furylhydrazone gives 2-methylene-2,3-dihydrofuran, a highly reactive ene in the ene reaction. J Org Chem. 2005 Apr 1;70(7):2862-5. PubMed PMID: 15787589.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481