


Cgs 15696

## Chemical Properties

CAS No. : 104615-00-1  
 Formula: C<sub>13</sub>H<sub>7</sub>ClN<sub>4</sub>O<sub>2</sub>  
 Molecular Weight: 286.67  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Cgs 15696 is an adenosine antagonist.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.4883 mL	17.4417 mL	34.8833 mL
5 mM	0.6977 mL	3.4883 mL	6.9767 mL
10 mM	0.3488 mL	1.7442 mL	3.4883 mL
50 mM	0.0698 mL	0.3488 mL	0.6977 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481