

Compound T129046

Chemical Properties


CAS No. :

Formula: C₂₀H₂₇N₃O₁₁

Molecular Weight: 457.43

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	2-phenyl-2-([3,4,5-trihydroxy-6-((3,4,5-trihydroxy-6-(hydroxymethyl)oxan-2-yl)oxy)methyl)oxan-2-yl]oxy)acetonitrile is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1861 mL	10.9306 mL	21.8613 mL
5 mM	0.4372 mL	2.1861 mL	4.3723 mL
10 mM	0.2186 mL	1.0931 mL	2.1861 mL
50 mM	0.0437 mL	0.2186 mL	0.4372 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481