


TMC 151C

Chemical Properties

CAS No. :	240807-75-4
Formula:	C ₄₁ H ₇₄ O ₁₅
Molecular Weight:	807.028
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	TMC 151C is a natural product that can be used as a reference standard. The CAS number of TMC 151C is 240807-75-4.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.2391 mL	6.1956 mL	12.3911 mL
5 mM	0.2478 mL	1.2391 mL	2.4782 mL
10 mM	0.1239 mL	0.6196 mL	1.2391 mL
50 mM	0.0248 mL	0.1239 mL	0.2478 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481