

(6E)-heptadeca-6,16-diene-1,2,4-triol

Chemical Properties

CAS No. :

Formula: C17H32O3

Molecular Weight: 284.44

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	(6E)-heptadeca-6,16-diene-1,2,4-triol is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.5157 mL	17.5784 mL	35.1568 mL
5 mM	0.7031 mL	3.5157 mL	7.0314 mL
10 mM	0.3516 mL	1.7578 mL	3.5157 mL
50 mM	0.0703 mL	0.3516 mL	0.7031 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481