


TC-AQP1-1

Chemical Properties

CAS No. : 37710-81-9
 Formula: C₁₂H₁₀O₄
 Molecular Weight: 218.21
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	TC-AQP1-1 is a blocker of the aquaporin 1 channel identified by virtual screening and inhibits water flux through hAQP1 in <i>Xenopus laevis</i> oocyte swelling assays at low micromolar concentrations.
Targets(IC50)	Aquaporin

Solubility Information

Solubility	DMSO: 50 mg/mL (229.14 mM), Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
In vivo Formulation	10% DMSO+90% Corn Oil: 2.5 mg/mL (11.46 mM), Sonication is recommended. <i>Please add the solvents sequentially, clarifying the solution as much as possible before adding the next one. Dissolve by heating and/or sonication if necessary. Working solution is recommended to be prepared and used immediately. The formulation provided above is for reference purposes only. In vivo formulations may vary and should be modified based on specific experimental conditions.</i>

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.5827 mL	22.9137 mL	45.8274 mL
5 mM	0.9165 mL	4.5827 mL	9.1655 mL
10 mM	0.4583 mL	2.2914 mL	4.5827 mL
50 mM	0.0917 mL	0.4583 mL	0.9165 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Seeliger D, et al. Discovery of novel human aquaporin-1 blockers. ACS Chem Biol. 2013 Jan 18;8(1):249-56.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481