

ETP-47037

Chemical Properties


CAS No. : 1316760-76-5

Formula: C₂₀H₂₇N₉O₃S

Molecular Weight: 473.55

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	ETP-47037 is a potent inhibitor of the PI3K α isoform (IC ₅₀ : 0.99 nM) and inhibits the PI3K β (IC ₅₀ : 49.2 nM), PI3K δ (IC ₅₀ : 7.13 nM) and PI3K γ (IC ₅₀ : 49.1 nM) isoforms. ETP-47037 has potential for telomere protection chemoregulation studies.
Targets(IC ₅₀)	Others,PI3K

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1117 mL	10.5585 mL	21.1171 mL
5 mM	0.4223 mL	2.1117 mL	4.2234 mL
10 mM	0.2112 mL	1.0559 mL	2.1117 mL
50 mM	0.0422 mL	0.2112 mL	0.4223 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481