

## Serratenediol diacetate

## Chemical Properties

CAS No. :	27832-84-4
Formula:	C <sub>34</sub> H <sub>54</sub> O <sub>4</sub>
Molecular Weight:	526.79
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

## Biological Description

Description	Serratenediol diacetate is a useful organic compound for research related to life sciences and the catalog number is T125132.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.8983 mL	9.4914 mL	18.9829 mL
5 mM	0.3797 mL	1.8983 mL	3.7966 mL
10 mM	0.1898 mL	0.9491 mL	1.8983 mL
50 mM	0.038 mL	0.1898 mL	0.3797 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481