

Compound N016-0050

Chemical Properties


CAS No. :

Formula: C73H90O21

Molecular Weight: 1303.5

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Compound N016-0050 is a useful organic compound for research related to life sciences and the catalog number is T131337.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.7672 mL	3.8358 mL	7.6717 mL
5 mM	0.1534 mL	0.7672 mL	1.5343 mL
10 mM	0.0767 mL	0.3836 mL	0.7672 mL
50 mM	0.0153 mL	0.0767 mL	0.1534 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481