


RU 24969 succinate

Chemical Properties

CAS No. : 107008-28-6
 Formula: C₁₈H₂₂N₂O₅
 Molecular Weight: 346.38
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	RU 24969 succinate, a 5-HT receptor agonist, exhibits K _i values of 0.38 nM for 5-HT 1B and 2.5 nM for 5-HT 1A, indicating its affinity towards these receptors. It reduces fluid consumption and enhances forward locomotion. This compound is utilized in neurological disease research.
Targets(IC50)	Others,5-HT Receptor

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.887 mL	14.435 mL	28.870 mL
5 mM	0.5774 mL	2.887 mL	5.774 mL
10 mM	0.2887 mL	1.4435 mL	2.887 mL
50 mM	0.0577 mL	0.2887 mL	0.5774 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481