

N-(Azido-PEG2)-N-Boc-PEG4-acid

Chemical Properties

CAS No. : 2093153-82-1


Formula: C₂₂H₄₂N₄O₁₀

Molecular Weight: 522.59

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	N-(Azido-PEG2)-N-Boc-PEG4-acid is a polyethylene glycol (PEG)-based linker developed for synthesizing PROTeolysis TArgeting Chimeras (PROTACs)[1].
Targets(IC50)	Others,PROTAC Linker
In vitro	PROTACs, composed of two ligands linked together—one targeting an E3 ubiquitin ligase and the other the target protein—utilize the intracellular ubiquitin-proteasome system to selectively degrade target proteins[1].

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9135 mL	9.5677 mL	19.1355 mL
5 mM	0.3827 mL	1.9135 mL	3.8271 mL
10 mM	0.1914 mL	0.9568 mL	1.9135 mL
50 mM	0.0383 mL	0.1914 mL	0.3827 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

An S, et al. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. EBioMedicine. 2018 Oct;36:553-562

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481