


Detirelix

Chemical Properties

CAS No. : 89662-30-6
 Formula: C₇₈H₁₀₅ClN₁₈O₁₃
 Molecular Weight: 1538.23
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	Detirelix is an LHRH antagonist.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.6501 mL	3.2505 mL	6.501 mL
5 mM	0.130 mL	0.6501 mL	1.3002 mL
10 mM	0.065 mL	0.325 mL	0.6501 mL
50 mM	0.013 mL	0.065 mL	0.130 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481