


2',3'-Epoxyindicolactone

Chemical Properties

| | |
|-------------------|---|
| CAS No. : | 312619-44-6 |
| Formula: | C ₂₁ H ₁₈ O ₇ |
| Molecular Weight: | 382.368 |
| Storage: | Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA. |



Biological Description

| | |
|---------------|---|
| Description | 2',3'-Epoxyindicolactone is a useful organic compound for research related to life sciences. The catalog number is T125072 and the CAS number is 312619-44-6. |
| Targets(IC50) | Others |

Preparing Stock Solutions

| | 1mg | 5mg | 10mg |
|-------|-----------|------------|------------|
| 1 mM | 2.6153 mL | 13.0763 mL | 26.1527 mL |
| 5 mM | 0.5231 mL | 2.6153 mL | 5.2305 mL |
| 10 mM | 0.2615 mL | 1.3076 mL | 2.6153 mL |
| 50 mM | 0.0523 mL | 0.2615 mL | 0.5231 mL |

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481