

Cottonseed oil

Chemical Properties

CAS No. : 8001-29-4

Formula:

Molecular Weight:

COTTONSEED OIL

Storage: Pure form: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Cottonseed oil, a cooking oil extracted from the seeds of cotton plants, has been considered the most insecticidal of vegetable oils.
Targets(IC50)	Others

Solubility Information

Solubility	DMSO: 0.9 mg/mL, Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481