

## Hydroxy-PEG5-Boc

## Chemical Properties

CAS No. : 850090-09-4


Formula: C17H34O8

Molecular Weight: 366.45

Keep away from direct sunlight

Storage: Pure form: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Hydroxy-PEG6-Boc, a PEG-based PROTAC linker, is utilized in PROTAC synthesis [1].
Targets(IC50)	PROTAC Linker
In vitro	PROTACs contain two distinct ligands connected by a linker: one ligand targets an E3 ubiquitin ligase, while the other targets a specific protein. By exploiting the intracellular ubiquitin-proteasome system, PROTACs selectively degrade target proteins [1].

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.7289 mL	13.6444 mL	27.2889 mL
5 mM	0.5458 mL	2.7289 mL	5.4578 mL
10 mM	0.2729 mL	1.3644 mL	2.7289 mL
50 mM	0.0546 mL	0.2729 mL	0.5458 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

## Reference

An S, et al. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. EBioMedicine. 2018 Oct;36:553-562

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481