

## Boc-PEG4-methyl propionate

## Chemical Properties

CAS No. : 2100306-74-7


Formula: C17H32O8

Molecular Weight: 364.435

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Boc-PEG4-methyl propionate is a polyethylene glycol (PEG)-based linker used in the synthesis of prodrug activation-targeting chimeras (PROTACs).
Targets(IC50)	Others,PROTAC Linker
In vitro	PROTACs consist of two distinct ligands connected by a linker: one ligand targets an E3 ubiquitin ligase, while the other targets the desired protein. They harness the intracellular ubiquitin-proteasome system to selectively degrade target proteins[1].

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.7439 mL	13.7197 mL	27.4394 mL
5 mM	0.5488 mL	2.7439 mL	5.4879 mL
10 mM	0.2744 mL	1.372 mL	2.7439 mL
50 mM	0.0549 mL	0.2744 mL	0.5488 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

## Reference

An S, et al. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. EBioMedicine. 2018 Oct;36:553-562

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481