


NLRP3 antagonist 1

Chemical Properties

CAS No. :	2454017-83-3
Formula:	C ₁₆ H ₁₈ N ₆ O
Molecular Weight:	310.35
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	NLRP3 antagonist 1 shows potential for the cancer research that is a potent NLRP3 antagonist. NLRP3 is involved in the body's intrinsic immunity against stress injury and pathogenic infections that is mainly expressed in macrophages and neutrophils [1].
Targets(IC50)	Others,NOD-like Receptor (NLR)

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.2222 mL	16.1108 mL	32.2217 mL
5 mM	0.6444 mL	3.2222 mL	6.4443 mL
10 mM	0.3222 mL	1.6111 mL	3.2222 mL
50 mM	0.0644 mL	0.3222 mL	0.6444 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481