


BGT1-IN-9

Chemical Properties

CAS No. : 185444-92-2
 Formula: C₅H₁₀ClN₃O₂
 Molecular Weight: 179.605
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

HCl



Biological Description

Description	BGT1-IN-9 is an M1 muscarinic agonist.
Targets(IC50)	Others,GABA Receptor

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.5676 mL	27.8381 mL	55.6762 mL
5 mM	1.1135 mL	5.5676 mL	11.1352 mL
10 mM	0.5568 mL	2.7838 mL	5.5676 mL
50 mM	0.1114 mL	0.5568 mL	1.1135 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Al-Khawaja A, Haugaard AS, Marek A, Löffler R, Thiesen L, Santiveri M, Damgaard M, Bundgaard C, Frølund B, Wellendorph P. Pharmacological Characterization of [(3)H]ATPCA as a Substrate for Studying the Functional Role of the Betaine/GABA Transporter 1 and the Creatine Transporter. *ACS Chem Neurosci*. 2018 Mar 21;9(3):545-554. doi: 10.1021/acscchemneuro.7b00351. Epub 2017 Nov 22. PubMed PMID: 29131576.

Al-Khawaja A, Petersen JG, Damgaard M, Jensen MH, Vogensen SB, Lie ME, Kragholm B, Bräuner-Osborne H, Clausen RP, Frølund B, Wellendorph P. Pharmacological identification of a guanidine-containing β -alanine analogue with low micromolar potency and selectivity for the betaine/GABA transporter 1 (BGT1). *Neurochem Res*. 2014 Oct;39(10):1988-96. doi: 10.1007/s11064-014-1336-9. Epub 2014 May 23. PubMed PMID: 24852577.

Suñol C, Babot Z, Cristòfol R, Sonnewald U, Waagepetersen HS, Schousboe A. A possible role of the non-GAT1 GABA transporters in transfer of GABA from GABAergic to glutamatergic neurons in mouse cerebellar neuronal cultures. *Neurochem Res*. 2010 Sep;35(9):1384-90. doi: 10.1007/s11064-010-0196-1. Epub 2010 May 30. PubMed PMID: 20512624.

Moeckel GW, Zhang L, Fogo AB, Hao CM, Pozzi A, Breyer MD. COX2 activity promotes organic osmolyte accumulation and adaptation of renal medullary interstitial cells to hypertonic stress. *J Biol Chem*. 2003 May 23;278(21):19352-7. Epub 2003 Mar 9. PubMed PMID: 12637551.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481