

Benzyloxy carbonyl-PEG3-NHS ester

Chemical Properties

CAS No. : 2100306-68-9


Formula: C₂₁H₂₇N₁O₉

Molecular Weight: 437.44

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Benzyloxy carbonyl-PEG3-NHS ester is a polyethylene glycol (PEG) derivative compound utilized as a PROTAC linker in the synthesis of PROTACs [1].
Targets(IC50)	Others,PROTAC Linker
In vitro	PROTACs consist of two ligands connected by a linker: one binds to an E3 ubiquitin ligase, and the other to the target protein. They utilize the intracellular ubiquitin-proteasome system to selectively degrade target proteins [1].

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.286 mL	11.4301 mL	22.8603 mL
5 mM	0.4572 mL	2.286 mL	4.5721 mL
10 mM	0.2286 mL	1.143 mL	2.286 mL
50 mM	0.0457 mL	0.2286 mL	0.4572 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

An S, et al. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. EBioMedicine. 2018 Oct;36:553-562

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481