

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1627 mL	10.8134 mL	21.6268 mL
5 mM	0.4325 mL	2.1627 mL	4.3254 mL
10 mM	0.2163 mL	1.0813 mL	2.1627 mL
50 mM	0.0433 mL	0.2163 mL	0.4325 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Chang JW, et al. Highly selective inhibitors of monoacylglycerol lipase bearing a reactive group that is bioisosteric with endocannabinoid substrates. Chem Biol. 2012 May 25;19(5):579-88.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481