

TCO-PEG3-amine

Chemical Properties

CAS No. : 1800507-93-0


Formula: C17H32N2O5

Molecular Weight: 344.452

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	TCO-PEG3-amine is a PEG-based linker for PROTACs that joins two essential ligands, crucial for forming PROTAC molecules. This linker enables selective protein degradation by leveraging the ubiquitin-proteasome system within cells.
Targets(IC50)	Others,PROTAC Linker
In vitro	PROTACs utilize two ligands linked together: one binds to an E3 ubiquitin ligase and the other to the target protein. By exploiting the intracellular ubiquitin-proteasome system, PROTACs selectively degrade target proteins.

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.9032 mL	14.5159 mL	29.0318 mL
5 mM	0.5806 mL	2.9032 mL	5.8064 mL
10 mM	0.2903 mL	1.4516 mL	2.9032 mL
50 mM	0.0581 mL	0.2903 mL	0.5806 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Wang Y, et al. Degradation of proteins by PROTACs and other strategies. Acta Pharm Sin B. 2020 Feb;10(2):207-238.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481