


## Silychristin B

## Chemical Properties

CAS No. :	879325-58-3
Formula:	C <sub>25</sub> H <sub>22</sub> O <sub>10</sub>
Molecular Weight:	482.44
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Silychristin B is a natural product for research related to life sciences. The catalog number is TN5023 and the CAS number is 879325-58-3.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0728 mL	10.364 mL	20.728 mL
5 mM	0.4146 mL	2.0728 mL	4.1456 mL
10 mM	0.2073 mL	1.0364 mL	2.0728 mL
50 mM	0.0415 mL	0.2073 mL	0.4146 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481