


KW 4099

## Chemical Properties

CAS No. : 141364-05-8  
 Formula: C<sub>29</sub>H<sub>30</sub>NNaO<sub>3</sub>S  
 Molecular Weight: 495.61  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	KW 4099 is a thromboxane A <sub>2</sub> antagonist.
Targets(IC <sub>50</sub> )	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0177 mL	10.0886 mL	20.1772 mL
5 mM	0.4035 mL	2.0177 mL	4.0354 mL
10 mM	0.2018 mL	1.0089 mL	2.0177 mL
50 mM	0.0404 mL	0.2018 mL	0.4035 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481