


## cIAP1 Ligand-Linker Conjugates 3

## Chemical Properties

CAS No. :	2222354-20-1
Formula:	C39H56N4O11S
Molecular Weight:	788.95
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



## Biological Description

Description	cIAP1 Ligand-Linker Conjugates 3 are chemical compounds comprising an IAP ligand designed specifically for the E3 ubiquitin ligase and a PROTAC linker. These conjugates, denoted as cIAP1 Ligand-Linker Conjugates 3, find utility in the design of SNIPERs, short for Specific and Non-genetic IAP-dependent Protein Erasers [1].
Targets(IC50)	Others,E3 Ligase Ligand-Linker Conjugates,PROTAC Linker

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.2675 mL	6.3375 mL	12.6751 mL
5 mM	0.2535 mL	1.2675 mL	2.535 mL
10 mM	0.1268 mL	0.6338 mL	1.2675 mL
50 mM	0.0254 mL	0.1268 mL	0.2535 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

## Reference

Heterocyclic compound. WO2018066545A1.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481