

MRTX-EX185

Chemical Properties

CAS No. :

Formula: C33H33FN6O2

Molecular Weight: 564.65

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	MRTX-EX185 is an effective inhibitor of GDP-loaded KRAS and KRAS (G12D) with an IC50=90 nM for KRAS (G12D). MRTX-EX185 is also capable of binding GDP-loaded HRAS.
Targets(IC50)	Others,Ras,Kras

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.771 mL	8.855 mL	17.7101 mL
5 mM	0.3542 mL	1.771 mL	3.542 mL
10 mM	0.1771 mL	0.8855 mL	1.771 mL
50 mM	0.0354 mL	0.1771 mL	0.3542 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481