

Meclizine

Chemical Properties


CAS No. : 569-65-3

Formula: C₂₅H₂₇ClN₂

Molecular Weight: 390.95

Storage: Pure form: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Meclizine, also known as Meclozine, is a piperazine class H1 antagonist with antihistamine properties that inhibits the interaction of histamine at H1 receptors reversibly. It serves as an effective anti-motion sickness agent, able to cross the blood-brain barrier. Additionally, Meclizine acts as an agonist ligand for the mouse constitutive androstane receptor (CAR) and an inverse agonist for human CAR, making it useful in research on polyQ toxicity disorders, including Huntington's disease.
Targets(IC50)	Apoptosis,Histamine Receptor
In vitro	Meclizine (Meclozine; 50 μM; 24 hours) notably enhances the survival of ST Hdh cells following serum deprivation for 24 hours, mainly through the inhibition of apoptosis, as indicated by the reduction in caspase 3 and 7 cleavage. This protective effect is dose-dependent, displaying an EC ₅₀ of 17.3 μm, and achieving a maximal increase in cell survival by 218% compared to the control. Additionally, Meclizine offers protection to both mutant (ST HdhQ111/111) and wild-type (ST HdhQ7/7) striatal cells harboring polyglutamine (polyQ)-expanded huntingtin proteins against apoptosis induced by serum withdrawal. This was demonstrated in studies employing murine striatal cells, where treatment with a 50 μM concentration for 24 hours significantly increased cell survival, and Western blot analysis confirmed apoptosis suppression through caspase 3 and 7 cleavage analysis.
In vivo	Meclizine (also known as Meclozine; administered at doses ranging from 10-100 mg/kg intraperitoneally) has been shown to protect mice against kidney ischemia. Specifically, a pretreatment regimen of 100 mg/kg of Meclizine 17 or 24 hours before induced ischemia notably enhances kidney protection in the subject mice. The protective mechanism of Meclizine involves the direct inhibition of the Kennedy pathway responsible for phosphatidylethanolamine biosynthesis, which leads to a reduction in mitochondrial oxygen consumption, coupled with an upregulation of glycolysis. This study was conducted on 8-10 week old male C57BL/6 mice, with various dosages tested, revealing that the compound effectively shields these mice from kidney ischemia-reperfusion injury.

Solubility Information

A DRUG SCREENING EXPERT

Solubility	DMSO: 4.5 mg/mL (11.51 mM), Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.5579 mL	12.7894 mL	25.5787 mL
5 mM	0.5116 mL	2.5579 mL	5.1157 mL
10 mM	0.2558 mL	1.2789 mL	2.5579 mL
50 mM	0.0512 mL	0.2558 mL	0.5116 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Priti N. Patel, et al. Safety and Efficacy in the Treatment and Prevention of Motion Sickness. *Clinical Medicine Insights: Therapeutics*. 2011;3.

Vishal M Gohil, et al. Meclizine is neuroprotective in models of Huntington's disease. *Hum Mol Genet*. 2011 Jan 15; 20(2):294-300.

Wendong Huang, et al. Meclizine Is an Agonist Ligand for Mouse Constitutive Androstane Receptor (CAR) and an Inverse Agonist for Human CAR. *Mol Endocrinol*. 2004 Oct;18(10):2402-8.

Seiji Kishi, et al. Meclizine Preconditioning Protects the Kidney Against Ischemia-Reperfusion Injury. *EBioMedicine*. 2015 Jul 29;2(9):1090-101.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481