


2-Deacetyltaxachitriene A

Chemical Properties

CAS No. :	214769-96-7
Formula:	C ₃₀ H ₄₂ O ₁₂
Molecular Weight:	594.654
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	2-Deacetyltaxachitriene A is a natural product for research related to life sciences. The catalog number is TN2753 and the CAS number is 214769-96-7.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.6817 mL	8.4083 mL	16.8166 mL
5 mM	0.3363 mL	1.6817 mL	3.3633 mL
10 mM	0.1682 mL	0.8408 mL	1.6817 mL
50 mM	0.0336 mL	0.1682 mL	0.3363 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481