

## Compound T128821

## Chemical Properties

CAS No. :

Formula: C<sub>21</sub>H<sub>36</sub>O<sub>11</sub>

Molecular Weight: 464.51

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	(2R,3R,4S,5S,6R)-2-(((2E)-6-hydroxy-2,6-dimethylocta-2,7-dien-1-yl]oxy)-6-(((2S,3R,4S,5S)-3,4,5-trihydroxyoxan-2-yl]oxy)methyl)oxane-3,4,5-triol is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1528 mL	10.764 mL	21.5281 mL
5 mM	0.4306 mL	2.1528 mL	4.3056 mL
10 mM	0.2153 mL	1.0764 mL	2.1528 mL
50 mM	0.0431 mL	0.2153 mL	0.4306 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481