

Compound T128280

Chemical Properties


CAS No. :

Formula: C19H26O12

Molecular Weight: 446.41

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	(2R,3R,4S,5S,6R)-2-([(2H-1,3-benzodioxol-5-yl)methoxy]-6-([(2R,3R,4R)-3,4-dihydroxy-4-(hydroxymethyl)oxolan-2-yl]oxy)methyl)oxane-3,4,5-triol is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2401 mL	11.2005 mL	22.4009 mL
5 mM	0.448 mL	2.2401 mL	4.4802 mL
10 mM	0.224 mL	1.120 mL	2.2401 mL
50 mM	0.0448 mL	0.224 mL	0.448 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481