


Beraprost

Chemical Properties

CAS No. :	88430-50-6
Formula:	C ₂₄ H ₃₀ O ₅
Molecular Weight:	398.49
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Beraprost is a stable prostacyclin analog.
Targets(IC50)	Others, Prostaglandin Receptor

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.5095 mL	12.5474 mL	25.0947 mL
5 mM	0.5019 mL	2.5095 mL	5.0189 mL
10 mM	0.2509 mL	1.2547 mL	2.5095 mL
50 mM	0.0502 mL	0.2509 mL	0.5019 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Zhang P, Xu L, Guan H, Liu L, Liu J, Huang Z, Cao X, Liao Z, Xiao H, Li Y. Beraprost sodium, a prostacyclin analogue, reduces fructose-induced hepatocellular steatosis in mice and in vitro via the microRNA-200a and SIRT1 signaling pathway. *Metabolism*. 2017 Aug;73:9-21. doi: 10.1016/j.metabol.2017.05.003. Epub 2017 May 11. PubMed PMID: 28732575.

Zhang H, Li X, Huang J, Li H, Su Z, Wang J. Comparative Efficacy and Safety of Prostacyclin Analogs for Pulmonary Arterial Hypertension: A Network Meta-Analysis. *Medicine (Baltimore)*. 2016 Jan;95(4):e2575. doi: 10.1097/MD.0000000000002575. Review. PubMed PMID: 26825901; PubMed Central PMCID: PMC5291571.

Del Pozo R, Hernandez Gonzalez I, Escribano-Subias P. The prostacyclin pathway in pulmonary arterial hypertension: a clinical review. *Expert Rev Respir Med*. 2017 Jun;11(6):491-503. doi: 10.1080/17476348.2017.1317599. Epub 2017 Apr 24. Review. PubMed PMID: 28399721.

Kim M, Kim JU, Kim SM, Kim H. Effectiveness of beraprost sodium in maintaining vascular access patency in patients on hemodialysis. *Int Urol Nephrol*. 2017 Jul;49(7):1287-1295. doi: 10.1007/s11255-017-1586-y. Epub 2017 Apr 13. PubMed PMID: 28409403.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481