


2-Chloroanthraquinone

Chemical Properties

CAS No. :	131-09-9
Formula:	C ₁₄ H ₇ ClO ₂
Molecular Weight:	242.66
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	2-Chloroanthraquinone is a biochemical.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.121 mL	20.605 mL	41.2099 mL
5 mM	0.8242 mL	4.121 mL	8.242 mL
10 mM	0.4121 mL	2.0605 mL	4.121 mL
50 mM	0.0824 mL	0.4121 mL	0.8242 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Kamijo S, Takao G, Kamijo K, Tsuno T, Ishiguro K, Murafuji T. Alkylation of Nonacidic C(sp³)-H Bonds by Photoinduced Catalytic Michael-Type Radical Addition. *Org Lett.* 2016 Oct 7;18(19):4912-4915. Epub 2016 Sep 26. PubMed PMID: 27668684.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481