

d-Laserpitin

Chemical Properties


CAS No. :

Formula: C₁₉H₂₀O₆

Molecular Weight: 344.363

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	d-Laserpitin is a useful organic compound for research related to life sciences and the catalog number is T126023.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.9039 mL	14.5197 mL	29.0394 mL
5 mM	0.5808 mL	2.9039 mL	5.8079 mL
10 mM	0.2904 mL	1.452 mL	2.9039 mL
50 mM	0.0581 mL	0.2904 mL	0.5808 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481