


2-Demethoxyleptostachyol acetate

Chemical Properties

CAS No. : 139405-55-3
 Formula: C₂₅H₂₆O₁₁
 Molecular Weight: 502.472
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

| | |
|-------------|--|
| Description | 2-Demethoxyleptostachyol acetate is a natural product for research related to life sciences. The catalog number is TN6167 and the CAS number is 139405-55-3. |
|-------------|--|

Preparing Stock Solutions

| | 1mg | 5mg | 10mg |
|-------|-----------|-----------|------------|
| 1 mM | 1.9902 mL | 9.9508 mL | 19.9017 mL |
| 5 mM | 0.398 mL | 1.9902 mL | 3.9803 mL |
| 10 mM | 0.199 mL | 0.9951 mL | 1.9902 mL |
| 50 mM | 0.0398 mL | 0.199 mL | 0.398 mL |

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481